

Directions by train

Epidauros Biotechnologie AG
Am Neuland 1
82347 Bernried
Phone: +49 (0)8158 9985-0
Internet: www.epidauros.com

From Munich Airport to Bernried via Munich Central Station ("München Hauptbahnhof")

Take "S-Bahn" S8 or S1 to Munich Central Station ("München Hauptbahnhof").
(duration appr. 45 minutes)
Change "S-Bahn" and take S6 (direction Tutzing) until Tutzing.
(duration appr. 45 minutes)
Change from S-Bahn to regular train (direction to Kochel). Next stop is Bernried.
(duration appr. 5 minutes)

Follow the footpath at the south end of the train platform according to the map shown below.

Direct regular trains from Munich may be faster although their frequency is lower.
For traveling schedules using Deutsche Bahn see the internet portal
www.bahn.com. => "Internat. Guests"

